

## Product Datasheet

### Human IL-1 beta ELISpot Kit 856.101.001PC

<b>Product Name</b>	Human IL-1 beta ELISpot Kit
<b>Catalog Number</b>	856.101.001PC
<b>Description</b>	Diaclone Pre-coated ELISpot kits include precoated PVDF plates, Detection antibody, Alkaline phosphatase conjugate, BSA, BCIP/NBT ready-to-use substrate buffer.
<b>Tested Applications</b>	ELISpot
<b>Alternative Names</b>	IL-1b, IL-1beta
<b>Brand</b>	Diaclone
<b>Storage</b>	Store at +4°C
<b>Cross Reactivity</b>	No cross reactivity with other human cytokines
<b>Specificity</b>	Recognizes natural (pro and mature) human IL-1b