

Product Datasheet

Human IL-1 beta ELISpot Kit 856.101.001PC

Product Name	Human IL-1 beta ELISpot Kit
Catalog Number	856.101.001PC
Description	Diaclone Pre-coated ELISpot kits include precoated PVDF plates, Detection antibody, Alkaline phosphatase conjugate, BSA, BCIP/NBT ready-to-use substrate buffer.
Tested Applications	ELISpot
Alternative Names	IL-1b, IL-1beta
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other human cytokines
Specificity	Recognizes natural (pro and mature) human IL-1b