

## Product Datasheet

### Human IL-1 beta ELISpot Set 856.101.005P

<b>Product Name</b>	Human IL-1 beta ELISpot Set
<b>Catalog Number</b>	856.101.005P
<b>Description</b>	Diaclone ELISpot Sets include capture and detection antibodies, Streptavidin - Alkaline Phosphatase conjugated, BSA, BCIP/NTB, blocking reagent, and 5 non-sterile plates.
<b>Tested Applications</b>	ELISpot
<b>Alternative Names</b>	IL-1b, IL-1beta
<b>Brand</b>	Diaclone
<b>Storage</b>	Store at +4°C
<b>Note</b>	Spots read visually or with reader
<b>Cross Reactivity</b>	No cross reactivity with other human cytokines
<b>Specificity</b>	Recognizes natural (pro and mature) human IL-1b