

Product Datasheet

Recombinant PCT antigen, 100 µg 610080

Product Name	Recombinant PCT antigen, 100 µg
Catalog Number	610080
Description	Recombinant human procalcitonin (PCT) protein with a C-terminal histidine tag.
Species/Host	E. coli
Brand	Medix Biochemica
Form/Appearance	Lyophilized
Storage	Storage at 2-8 °C. After reconstitution stable for 1 week at 2-8 °C and 10 months at - 20 °C.
Associated Products	Anti-h PCT 4003 Anti-h PCT 4004 Anti-h PCT 4005 Anti-h PCT 4006 Anti-h PCT 4008 Anti-h Calcitonin 4051
Buffer	Lyophilized from 50 mM Tris-HCl, pH 7.5; 150 mM NaCl; containing 6 % sucrose as stabilizer
Molecular Weight	Predicted molecular weight: 14 kDa
Shelf Life	Unspecified for lyophilized product