

## Product Datasheet

### Human IL-5 ELISpot Set 856.111.001P

<b>Product Name</b>	Human IL-5 ELISpot Set
<b>Catalog Number</b>	856.111.001P
<b>Description</b>	Diaclone ELISpot Sets include capture and detection antibodies, Streptavidin - Alkaline Phosphatase conjugated, BSA, BCIP/NTB, blocking reagent, and 1 non-sterile plate.
<b>Tested Applications</b>	ELISpot
<b>Alternative Names</b>	NULL
<b>Brand</b>	Diaclone
<b>Storage</b>	Store at +4°C
<b>Note</b>	Spots read visually or with reader
<b>Cross Reactivity</b>	No cross reactivity with other human cytokines
<b>Specificity</b>	Recognizes both natural and recombinant human IL-5