

Product Datasheet

Recombinant HRP2 antigen, 100 µg 610081

Product Name	Recombinant HRP2 antigen, 100 µg
Catalog Number	610081
Description	Recombinant Plasmodium falciparum histidine-rich protein 2 (HRP2).
Species/Host	E. coli
Brand	Medix Biochemica
Form/Appearance	Lyophilized
Storage	Storage at 2-8 °C. After reconstitution stable for 1 month at 2-8 °C and 1 months at - 20 °C.
Associated Products	Anti-HRP2 3901, Anti-HRP2 3902
Buffer	Lyophilized from 25 mM sodium acetate, pH 5.0, 0.02 % Tween 20, 150 mM NaCl; containing 6 % sucrose as a stabilizer
Molecular Weight	Predicted molecular weight: 33 kDa
Shelf Life	Unspecified for lyophilized product