

Product Datasheet

Recombinant NT-proBNP antigen, 100 µg 610090

Product Name	Recombinant NT-proBNP antigen, 100 µg
Catalog Number	610090
Description	Recombinant human N-terminal pro-B-type natriuretic peptide (NT-proBNP) with an N-terminal histidine tag.
Species/Host	E. coli
Brand	Medix Biochemica
Form/Appearance	Lyophilized
Storage	Storage at 2-8°C. After reconstitution stable for 1 month at 2-8°C and 1 month at -20°C.
Associated Products	Anti-h NT-proBNP 1306/RC1306 Anti-h NT-proBNP 1307/RC1307 Anti-h NT-proBNP 1308 Anti-h NT-proBNP 1309 Anti-h NT-proBNP 1310 Anti-h NT-proBNP 1311 Anti-h NT-proBNP 1312
Buffer	Lyophilized from 50 mM Tris-HCl, pH 7.5; 150 mM NaCl; containing 3 % sucrose, 2 % D-mannitol and 0.01 % Tween 20 as stabilizers.
Molecular Weight	Predicted molecular weight: 11 kDa
Shelf Life	Unspecified for lyophilized product