

## Product Datasheet

### Recombinant NT-proBNP antigen, 100 µg 610090

<b>Product Name</b>	Recombinant NT-proBNP antigen, 100 µg
<b>Catalog Number</b>	610090
<b>Description</b>	Recombinant human N-terminal pro-B-type natriuretic peptide (NT-proBNP) with an N-terminal histidine tag.
<b>Species/Host</b>	E. coli
<b>Brand</b>	Medix Biochemica
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	Storage at 2-8°C. After reconstitution stable for 1 month at 2-8°C and 1 month at -20°C.
<b>Associated Products</b>	Anti-h NT-proBNP 1306/RC1306 Anti-h NT-proBNP 1307/RC1307 Anti-h NT-proBNP 1308 Anti-h NT-proBNP 1309 Anti-h NT-proBNP 1310 Anti-h NT-proBNP 1311 Anti-h NT-proBNP 1312
<b>Buffer</b>	Lyophilized from 50 mM Tris-HCl, pH 7.5; 150 mM NaCl; containing 3 % sucrose, 2 % D-mannitol and 0.01 % Tween 20 as stabilizers.
<b>Molecular Weight</b>	Predicted molecular weight: 11 kDa
<b>Shelf Life</b>	Unspecified for lyophilized product