

Product Datasheet

Human IL-13 ELISpot Set 856.121.001P

Product Name	Human IL-13 ELISpot Set
Catalog Number	856.121.001P
Description	Diaclone ELISpot Sets include capture and detection antibodies, Streptavidin - Alkaline Phosphatase conjugated, BSA, BCIP/NTB, blocking reagent, and 1 non-sterile plate.
Tested Applications	ELISpot
Alternative Names	NULL
Brand	Diaclone
Storage	Store at +4°C
Note	Spots read visually or with reader
Cross Reactivity	No cross reactivity with other human cytokines
Specificity	Recognizes both natural and recombinant human IL-13