

Product Datasheet

Recombinant Cystatin C antigen, 100 µg 610100

Product Name	Recombinant Cystatin C antigen, 100 µg
Catalog Number	610100
Description	Recombinant human cystatin C antigen with an N-terminal histidine tag.
Species/Host	E. coli
Brand	Medix Biochemica
Form/Appearance	Lyophilized
Storage	Storage at 2-8 °C. After reconstitution 1 month at 2-8 °C and 1 month at - 20 °C.
Associated Products	Anti-h Cystatin C 10001, 10002, 10004 and 10005
Buffer	Lyophilized from 50 mM Tris-HCl, pH 7.5, 150 mM NaCl; containing 6 % sucrose as a stabilizer
Molecular Weight	Predicted molecular weight: 16 kDa
Shelf Life	Unspecified for lyophilized product