

Product Datasheet

Recombinant cTnT antigen, 100 µg 610101

Product Name	Recombinant cTnT antigen, 100 µg
Catalog Number	610101
Description	Recombinant human cardiac troponin T antigen (isoform 6) with a C-terminal histidine tag.
Species/Host	E. coli
Alternative Names	Troponin T
Brand	Medix Biochemica
Form/Appearance	Lyophilized
Storage	Storage at 2-8 °C. After reconstitution stable for 1 month at 2-8 °C and 6 months at - 20 °C.
Associated Products	Anti-h cTnT 3701 Anti-h cTnT 3703 Anti-h cTnT 3708 Anti-h cTnT 3710 Anti-h cTnT 3711 Anti-h cTnT 3712
Buffer	Lyophilized from 50 mM Tris-HCl, pH 7.5; 150 mM NaCl, 4 M urea; containing 3 % sucrose, 2 % D-mannitol, and 0.01 % Tween 20 as stabilizers
Molecular Weight	Predicted molecular weight: 36 kDa
Shelf Life	Unspecified for lyophilized product