

Product Datasheet

Recombinant cTnI antigen, 100 µg 610102

Product Name	Recombinant cTnI antigen, 100 µg
Catalog Number	610102
Description	Recombinant human cardiac troponin I antigen with a C-terminal histidine tag.
Species/Host	E. coli
Alternative Names	troponin
Brand	Medix Biochemica
Form/Appearance	Lyophilized
Storage	+2-8°C
Associated Products	Anti-h cTnI 9701/RC970, Anti-h cTnI 9703, Anti-h cTnI 9705, Anti-h cTnI 9707/RC9707, Anti-h cTnI RC9750
Buffer	Lyophilized from 50 mM Tris-HCl, pH 7.5, 150 mM NaCl, 8 M urea; containing 3 % sucrose, 2 % D-mannitol, and 0.01 % Tween 20 as stabilizers
Molecular Weight	Predicted molecular weight: 26 kDa
Shelf Life	Unspecified