

Product Datasheet

Human Perforin ELISpot Kit 856.141.005PC

| | |
|----------------------------|--|
| Product Name | Human Perforin ELISpot Kit |
| Catalog Number | 856.141.005PC |
| Description | Diaclone Pre-coated ELISpot kits include precoated PVDF plates, Detection antibody, Alkaline phosphatase conjugate, BSA, BCIP/NBT ready-to-use substrate buffer. |
| Tested Applications | ELISpot |
| Alternative Names | NULL |
| Brand | Diaclone |
| Storage | Store at +4°C |
| Specificity | Recognizes natural human Perforin |