

Product Datasheet

Human IL-17F ELISpot Kit 856.161.001PC

| | |
|----------------------------|--|
| Product Name | Human IL-17F ELISpot Kit |
| Catalog Number | 856.161.001PC |
| Description | Diaclone Pre-coated ELISpot kits include precoated PVDF plates, Detection antibody, Alkaline phosphatase conjugate, BSA, BCIP/NBT ready-to-use substrate buffer. |
| Tested Applications | ELISpot |
| Alternative Names | NULL |
| Brand | Diaclone |
| Storage | Store at +4°C |
| Specificity | Recognizes natural human IL-17F |