

Product Datasheet

Recombinant Galectin-3 antigen, 100 µg 610144

| | |
|----------------------------|---|
| Product Name | Recombinant Galectin-3 antigen, 100 µg |
| Catalog Number | 610144 |
| Description | Recombinant human Galectin-3 antigen with an N-terminal histidine tag. |
| Species/Host | E. coli |
| Brand | Medix Biochemica |
| Form/Appearance | Liquid |
| Concentration | Lot specific |
| Storage | +2-8°C, do not freeze |
| Associated Products | Anti-h Galectin-3 10301 Anti-h Galectin-3 10302 Anti-h Galectin-3 10303 Anti-h Galectin-3 10304 Anti-h Galectin-3 10305 |
| Buffer | 50 mM Tris-HCl, 150 mM NaCl, 1 mM TCEP, pH 7.5, 0.095 % NaN ₃ as a preservative |
| Molecular Weight | Predicted molecular weight: 28 kDa |
| Shelf Life | Unspecified |