

Product Datasheet

Human IL-17A/F ELISpot Kit 856.171.002PC

Product Name	Human IL-17A/F ELISpot Kit
Catalog Number	856.171.002PC
Description	Diaclone Pre-coated ELISpot kits include precoated PVDF plates, Detection antibody, Alkaline phosphatase conjugate, BSA, BCIP/NBT ready-to-use substrate buffer.
Tested Applications	ELISpot
Alternative Names	NULL
Brand	Diaclone
Storage	Store at +4°C
Specificity	Recognizes natural human IL-17A/F