

## Product Datasheet

### Recombinant NSE antigen, 100 µg 610150

<b>Product Name</b>	Recombinant NSE antigen, 100 µg
<b>Catalog Number</b>	610150
<b>Description</b>	Recombinant human neuron-specific enolase (NSE) γ-isoform with a C-terminal histidine tag.
<b>Species/Host</b>	E. coli
<b>Brand</b>	Medix Biochemica
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	Lot specific
<b>Storage</b>	+2-8°C
<b>Associated Products</b>	<b>Anti-h NSE 9601</b> <b>Anti-h NSE 9602</b>
<b>Buffer</b>	50 mM Tris-HCl, 150 mM NaCl, 5 mM MgSO <sub>4</sub> , 0.02 % Tween 20, pH 7.5, 0.095 % NaN <sub>3</sub> as a preservative
<b>Molecular Weight</b>	Predicted molecular weight: 48 kDa
<b>Shelf Life</b>	Unspecified