

Product Datasheet

Recombinant NSE antigen, 100 µg 610150

Product Name	Recombinant NSE antigen, 100 µg
Catalog Number	610150
Description	Recombinant human neuron-specific enolase (NSE) γ-isoform with a C-terminal histidine tag.
Species/Host	E. coli
Brand	Medix Biochemica
Form/Appearance	Liquid
Concentration	Lot specific
Storage	+2-8°C
Associated Products	Anti-h NSE 9601 Anti-h NSE 9602
Buffer	50 mM Tris-HCl, 150 mM NaCl, 5 mM MgSO ₄ , 0.02 % Tween 20, pH 7.5, 0.095 % NaN ₃ as a preservative
Molecular Weight	Predicted molecular weight: 48 kDa
Shelf Life	Unspecified