

## Product Datasheet

### Murine IL-2 ELISpot Set 862.001.005S

<b>Product Name</b>	Murine IL-2 ELISpot Set
<b>Catalog Number</b>	862.001.005S
<b>Description</b>	Diaclone ELISpot Sets include capture and detection antibodies, Streptavidin - Alkaline Phosphatase conjugated, BSA, BCIP/NTB, blocking reagent, and 5 sterile plates.
<b>Tested Applications</b>	ELISpot
<b>Alternative Names</b>	NULL
<b>Brand</b>	Diaclone
<b>Storage</b>	Store at +4°C
<b>Note</b>	Spots read visually or with reader
<b>Cross Reactivity</b>	No cross reactivity with other murine cytokines
<b>Specificity</b>	Recognizes natural murine IL-2