

Product Datasheet

Murine IFN gamma ELISpot Kit 862.031.002PC

Product Name	Murine IFN gamma ELISpot Kit
Catalog Number	862.031.002PC
Description	Diaclone Pre-coated ELISpot kits include precoated PVDF plates, Detection antibody, Alkaline phosphatase conjugate, BSA, BCIP/NBT ready-to-use substrate buffer.
Tested Applications	ELISpot
Alternative Names	IFN-g, IFN-gamma
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other murine cytokines
Specificity	Recognizes natural murine IFN-g