

Product Datasheet

Murine IFN gamma ELISpot Set 862.031.005P

Product Name	Murine IFN gamma ELISpot Set
Catalog Number	862.031.005P
Description	Diaclone ELISpot Sets include capture and detection antibodies, Streptavidin - Alkaline Phosphatase conjugated, BSA, BCIP/NTB, blocking reagent, and 5 non-sterile plates.
Tested Applications	ELISpot
Alternative Names	IFN-g, IFN-gamma
Brand	Diaclone
Storage	Store at +4°C
Note	Spots read visually or with reader
Cross Reactivity	No cross reactivity with other murine cytokines
Specificity	Recognizes natural murine IFN-g