

Product Datasheet

Anti-h Insulin 9350 700013

Product Name	Anti-h Insulin 9350
Catalog Number	700013
Description	Monoclonal antibody raised against human Insulin. In vitro cultured mouse monoclonal antibody.
Tested Applications	ELISA
Brand	Medix Biochemica
Form/Appearance	Liquid
Concentration	> 1.0 mg/ml
Storage	+2-8°C
Isotype	IgG1
Clonality	Monoclonal
Epitope	N/D
Purity	≥ 95 %
Affinity constant	N/D
Buffer	0.1 M PBS, pH 7.4, 0.9 % NaCl, 0.05 % NaN ₃ as a preservative
IEF Profile	6.8 - 7.2
Cross Reactivity	N/D
Specificity	Antibody recognizes human insulin
Matched Pairs	Anti-h Insulin 9351 (700014)
Shelf Life	10 years