

## Product Datasheet

### Rat IFN gamma ELISpot Set 867.011.001P

<b>Product Name</b>	Rat IFN gamma ELISpot Set
<b>Catalog Number</b>	867.011.001P
<b>Description</b>	Diaclone ELISpot Sets include capture and detection antibodies, Streptavidin - Alkaline Phosphatase conjugated, BSA, BCIP/NTB, blocking reagent, and 1 non-sterile plate.
<b>Tested Applications</b>	ELISpot
<b>Alternative Names</b>	IFN-g, IFN-gamma
<b>Brand</b>	Diaclone
<b>Storage</b>	Store at +4°C
<b>Note</b>	Spots read visually or with reader
<b>Cross Reactivity</b>	No cross reactivity with other rat cytokines
<b>Specificity</b>	Recognizes natural rat IFN-g