

Product Datasheet

Rat IFN gamma ELISpot Set 867.011.001S

Product Name	Rat IFN gamma ELISpot Set
Catalog Number	867.011.001S
Description	Diaclone ELISpot Sets include capture and detection antibodies, Streptavidin - Alkaline Phosphatase conjugated, BSA, BCIP/NTB, blocking reagent, and 1 sterile plate.
Tested Applications	ELISpot
Alternative Names	IFN-g, IFN-gamma
Brand	Diaclone
Storage	Store at +4°C
Note	Spots read visually or with reader
Cross Reactivity	No cross reactivity with other rat cytokines
Specificity	Recognizes natural rat IFN-g