

Product Datasheet

Rat IFN gamma ELISpot Kit 867.011.005PC

Product Name	Rat IFN gamma ELISpot Kit
Catalog Number	867.011.005PC
Description	Diaclone Pre-coated ELISpot kits include precoated PVDF plates, Detection antibody, Alkaline phosphatase conjugate, BSA, BCIP/NBT ready-to-use substrate buffer.
Tested Applications	ELISpot
Alternative Names	IFN-g, IFN-gamma
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other rat cytokines
Specificity	Recognizes natural rat IFN-g