

Product Datasheet

Anti-h T3 10550 700015

Product Name	Anti-h T3 10550
Catalog Number	700015
Description	Monoclonal antibody raised against human T3. In vitro cultured mouse hybridomas.
Tested Applications	ELISA
Brand	Medix Biochemica
Form/Appearance	Liquid
Concentration	> 1.0 mg/ml
Storage	+2-8°C
Isotype	IgG1
Clonality	Monoclonal
Epitope	N/D
Purity	≥ 95 %
Affinity constant	N/D
Associated Products	Triiodothyronine (T3) antigen - BSA Conjugate 581-05
Buffer	0.1 M PBS, pH 7.4, 0.9 % NaCl, 0.05 % NaN ₃ as a preservative
IEF Profile	6.7-7.1
Cross Reactivity	L-3,3',5-triiodothyronine (L-T3) 100 % L-thyroxine (L-T4) 0,17 % D-thyroxine (D-T4) 0,04 % 3,3',5-triiodothyroacetic acid (TRIAC) 52 % 3,5-diiodo-L-tyrosine 0,22 %
Specificity	Antibody recognizes human triiodothyronine (T3)
Shelf Life	10 years