

Product Datasheet

Anti-Vitamin 25 OH D2/D3 3850 700016

Product Name	Anti-Vitamin 25 OH D2/D3 3850
Catalog Number	700016
Description	Monoclonal antibody recognizing both 25 OH Vitamin D2 and 25 OH D3. In vitro cultured mouse monoclonal antibody.
Tested Applications	ELISA
Alternative Names	Vitamin D
Brand	Medix Biochemica
Form/Appearance	Liquid
Concentration	> 1.0 mg/ml
Storage	+2-8°C
Isotype	IgG1
Clonality	Monoclonal
Epitope	N/D
Affinity constant	N/D
Associated Products	Vitamin D antigen, BSA conjugated 710048 Vitamin D antigen, Biotin conjugated 710049 Vitamin D antigen, NH2 710070 Vitamin D antigen, COOH 710071
Buffer	0.1 M PBS, pH 7.4, 0.9 % NaCl, 0.05 % NaN3 as a preservative
IEF Profile	N/D
Cross Reactivity	Vitamin D3 < 3 %, Vitamin D2 ~ 1 %, 1, 25(OH)2 Vitamin D3* 20 %, 1, 25(OH)2 Vitamin D2* ~ 2 % . *The concentration in blood is very low, about 0.1 %
Specificity	Antibody recognizes natural 25 OH Vitamin D2 and 25 OH Vitamin D3
Shelf Life	10 years