

Product Datasheet

Anti-h Calcitonin 4051 700020

Product Name	Anti-h Calcitonin 4051
Catalog Number	700020
Description	Monoclonal antibody raised against human calcitonin. In vitro cultured mouse hybridomas.
Tested Applications	ELISA
Brand	Medix Biochemica
Form/Appearance	Liquid
Concentration	> 1.0 mg/ml
Storage	+2-8°C
Isotype	IgG1
Clonality	Monoclonal
Purity	≥ 95 %
Affinity constant	N/D
Associated Products	Recombinant PCT antigen, Medix Biochemica 610080
Buffer	0.1 M BBS, pH 8, 0.9 % NaCl
IEF Profile	6.0 - 6.3
Specificity	Antibody recognizes human calcitonin and procalcitonin
Shelf Life	10 years