

---

## Product Datasheet

### Human IFN $\gamma$ / IL-5 Dual ELISpot Set 874.020.005P

|                            |   |
|----------------------------|---|
| <b>Product Name</b>        | Human IFN $\gamma$ / IL-5 Dual ELISpot Set  |
| <b>Catalog Number</b>      | 874.020.005P  |
| <b>Description</b>         | Diaclone Dual ELISpot Sets include capture and detection antibodies for both of the targeted cytokines, Streptavidin-Alkaline Phosphatase conjugated, anti-FITC HRP conjugated mAb, BSA, BCIP/NTB and AEC, Blocking reagent and 5 non-sterile plates. |
| <b>Tested Applications</b> | ELISpot   |
| <b>Brand</b>               | Diaclone  |
| <b>Storage</b>             | Store at +4°C   |
| <b>Cross Reactivity</b>    | Cross-reactive with simian IFN-g  |
| <b>Specificity</b>         | Recognizes both natural and recombinant human IFN-g and IL-5  |