

Product Datasheet

Anti-h Renin 10850 700025

Product Name	Anti-h Renin 10850
Catalog Number	700025
Description	Monoclonal antibody raised against human renin. In vitro cultured mouse hybridomas.
Brand	Medix Biochemica
Form/Appearance	Liquid
Concentration	> 1.0 mg/ml
Storage	+2-8°C
Isotype	IgG1
Clonality	Monoclonal
Epitope	N/D
Buffer	0.1 M PBS, pH 7.4, 0.9 % NaCl, 0.05 % NaN ₃ as a preservative
IEF Profile	N/D
Cross Reactivity	Prorenin < 10 %
Specificity	Antibody recognizes human renin
Shelf Life	10 years