

Product Datasheet

Human IFN γ / Granzyme B Dual ELISpot Set 874.060.001P

Product Name	Human IFN γ / Granzyme B Dual ELISpot Set
Catalog Number	874.060.001P
Description	Diaclone Dual ELISpot Sets include capture and detection antibodies for both of the targeted cytokines, Streptavidin-Alkaline Phosphatase conjugated, anti-FITC HRP conjugated mAb, BSA, BCIP/NTB and AEC, Blocking reagent and 1 non-sterile plate.
Tested Applications	ELISpot
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other human cytokines. Cross reactivity with simian IFN-g
Specificity	Recognizes natural human IFN-g and Granzyme B