

Product Datasheet

Human IFN γ / IL-17A Dual ELISpot Set 874.080.001S

Product Name	Human IFN γ / IL-17A Dual ELISpot Set
Catalog Number	874.080.001S
Description	Diaclone Dual ELISpot Sets include capture and detection antibodies for both of the targeted cytokines, Streptavidin-Alkaline Phosphatase conjugated, anti-FITC HRP conjugated mAb, BSA, BCIP/NTB and AEC, blocking reagent and 1 sterile plate.
Tested Applications	ELISpot
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other human cytokines. Cross reactivity with simian IFN-g.
Specificity	Recognizes natural human IL-17A and IFN-g