

Product Datasheet

Human IFNγ / IL-17A Dual Fluorospot Set 874.082.005P

Product Name	Human IFNγ / IL-17A Dual Fluorospot Set
Catalog Number	874.082.005P
Description	Diaclone Dual Fluorospot Sets include capture antibody for cytokine 1, FITC-conjugated detection antibody for cytokine 1, anti-FITC antibody green fluorescence conjugate, capture antibody for cytokine 2, biotinylated detection antibody for cytokine 2, streptavidin-phycoerythrin conjugate, BSA. Kit can be sold with or without plates.
Tested Applications	Fluorospot
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other human cytokines
Specificity	Recognizes natural human IFN-γ and IL-17A