

Product Datasheet

Vitamin D antigen, BSA conjugated, 1 mg 710048

Product Name	Vitamin D antigen, BSA conjugated, 1 mg
Catalog Number	710048
Description	Vitamin D antigen, BSA conjugated, 1 mg
Brand	Medix Biochemica
Form/Appearance	Liquid, turbid. Vortex before use.
Concentration	> 0.5 mg/ml
Storage	Short time storage at -20 °C, long time storage at -80 °C.
Source	Synthetic
Associated Products	Anti-Vitamin 25 OH D2/D3 3850 Anti-Vitamin 25 OH D2/D3 3851
Buffer	0.9 % NaCl
Shelf Life	60 months