

Product Datasheet

Vitamin D antigen, Biotin conjug, 1 mg 710049

Product Name	Vitamin D antigen, Biotin conjug, 1 mg
Catalog Number	710049
Description	Vitamin D antigen, Biotin conjug, 1 mg
Brand	Medix Biochemica
Form/Appearance	Lyophilized antigen, on the bottom of the vial. Not visually detectable.
Concentration	1 mg/ml when reconstituted with 1000 µl of EtOH
Storage	Storage at -20 °C. After reconstitution 6 months at -20 °C. Long term storage at -80 °C.
Source	Synthetic
Associated Products	Anti-Vitamin 25 OH D2/D3 3850 Anti-Vitamin 25 OH D2/D3 3851
Shelf Life	60 months