

## Product Datasheet

### E2 (Estradiol) 6-HS antigen, 1 mg 710050

<b>Product Name</b>	E2 (Estradiol) 6-HS antigen, 1 mg
<b>Catalog Number</b>	710050
<b>Description</b>	Synthetic Estradiol 6-hemisuccinate antigen, COOH functionalized
<b>Brand</b>	Medix Biochemica
<b>Form/Appearance</b>	Lyophilized
<b>Concentration</b>	1 mg/ml when reconstituted with 1000 µl of acetone
<b>Storage</b>	Storage at -20 °C. After reconstitution stable for 6 months, storage at -20 °C.
<b>Associated Products</b>	Anti-h E2 (Estradiol) 10450
<b>Shelf Life</b>	60 months