

Product Datasheet

Vitamin D antigen, COOH, 1 mg 710071

Product Name Vitamin D antigen, COOH, 1 mg

Catalog Number 710071

Description Vitamin D antigen, COOH, 1 mg

Brand Medix Biochemica

Form/Appearance Lyophilized

Concentration 1 mg/ml when reconstituted with 1000 µl of EtOH

Storage Storage at -20 °C. After reconstitution 6 months at -20 °C. Long term storage at -80 °C.

Source Synthetic

Associated Products Anti-Vitamin 25 OH D2/D3 3850 Anti-Vitamin 25 OH D2/D3 3851

Shelf Life 60 months