

Product Datasheet

Anti-h Albumin 6502 SP-2 100064

Product Name	Anti-h Albumin 6502 SP-2
Catalog Number	100064
Description	Monoclonal mouse antibody, cultured in vitro under conditions free from animal-derived components.
Tested Applications	LF, FIA
Brand	Medix Biochemica
Form/Appearance	Liquid, may turn slightly opaque during storage
Concentration	2.5 mg/ml (+/- 10 %)
Storage	+2-8°C
Isotype	IgG2a
Clonality	Monoclonal
Epitope	N/D
Purity	≥ 95 %
Affinity constant	KA= 5 x 10 ⁹ 1/M
Buffer	0.9 % NaCl, 0.095 % NaN ₃ as a preservative
IEF Profile	6.9-9.2
Cross Reactivity	AFP < 0.7 %
Specificity	Antibody recognizes human albumin
Shelf Life	36 months