

Product Datasheet

Recombinant SARS-CoV-2 S1, 1 mg 710091

Product Name	Recombinant SARS-CoV-2 S1, 1 mg
Catalog Number	710091
Description	Recombinant human coronavirus SARS-CoV-2 Spike S1 protein contains amino acids 14-681 of the S protein, and two extra amino acids (AS) at the N-terminus and a His-6 tag at the C-terminus.
Species/Host	CHO Cells
Brand	Medix Biochemica
Form/Appearance	Liquid
Concentration	0.5-1.0 mg/ml
Storage	-70 °C
Purity	SDS-PAGE > 90 %
Buffer	PBS, pH 7.4
Molecular Weight	Predicted molecular weight: 76.4 kDa
Shelf Life	Unspecified