

## Product Datasheet

### Recombinant SARS-CoV-2 S1, 1 mg 710091

<b>Product Name</b>	Recombinant SARS-CoV-2 S1, 1 mg
<b>Catalog Number</b>	710091
<b>Description</b>	Recombinant human coronavirus SARS-CoV-2 Spike S1 protein contains amino acids 14-681 of the S protein, and two extra amino acids (AS) at the N-terminus and a His-6 tag at the C-terminus.
<b>Species/Host</b>	CHO Cells
<b>Brand</b>	Medix Biochemica
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	0.5-1.0 mg/ml
<b>Storage</b>	-70 °C
<b>Purity</b>	SDS-PAGE > 90 %
<b>Buffer</b>	PBS, pH 7.4
<b>Molecular Weight</b>	Predicted molecular weight: 76.4 kDa
<b>Shelf Life</b>	Unspecified