

Product Datasheet

Anti-Mouse IL-23 Monoclonal Antibody, Azide Free Clone B-B23 400.000

Product Name	Anti-Mouse IL-23 Monoclonal Antibody, Azide Free Clone B-B23
Catalog Number	400.000
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	ELISA
Immunogen	Recombinant murine IL-23p19
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Specificity	Recognises both natural and recombinant protein
Hybridoma	Myeloma X63/AG.8653 x C57 Black/6 spleen cells