

Product Datasheet

Human TGF beta 1 ELISA Kit 650.010

| | |
|----------------------------|--|
| Product Name | Human TGF beta 1 ELISA Kit |
| Catalog Number | 650.010 |
| Description | Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure. |
| Range | 31.3-2000 pg/ml |
| Sensitivity | 9 pg/ml |
| Tested Applications | ELISA |
| Alternative Names | TGFb1, TGF-beta1 |
| Brand | Diaclone |
| Storage | Store at +4°C |
| Cross Reactivity | Cross reactivity with sheep and pig TGF-β1 |
| Specificity | Recognizes both natural and recombinant human TGF-β1 |
| Sample Types | Serum, Plasma, Cell culture supernatant |