

Product Datasheet

Human VEGF R1 ELISA Kit 650.030

Product Name	Human VEGF R1 ELISA Kit
Catalog Number	650.030
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	0.16-10 ng/ml
Sensitivity	0.1 ng/ml
Tested Applications	ELISA
Brand	Diaclone
Storage	Store at +4°C
Specificity	Recognizes soluble VEGF-R1
Sample Types	Serum, Plasma, Cell culture supernatant