

Product Datasheet

Human Cytochrome C ELISA Kit 650.070

Product Name	Human Cytochrome C ELISA Kit
Catalog Number	650.070
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	0.16-5 ng/ml
Sensitivity	0.08 ng/ml
Tested Applications	ELISA
Alternative Names	CYC
Brand	Diaclone
Storage	Store at +4°C
Specificity	Specific for human Cytochrome C
Sample Types	Whole blood, Cell lysates