

Product Datasheet

Anti-h FSH 6602 SP-5 100067

Product Name	Anti-h FSH 6602 SP-5
Catalog Number	100067
Description	Monoclonal mouse antibody, cultured in vitro under conditions free from animal-derived components.
Tested Applications	FIA
Brand	Medix Biochemica
Form/Appearance	Liquid, may turn slightly opaque during storage
Concentration	5.0 mg/ml (+/- 10 %)
Storage	+2-8°C
Isotype	IgG1
Clonality	Monoclonal
Epitope	N/D
Purity	≥ 95 %
Affinity constant	3×10^{10} 1/M
Associated Products	Antibody recognizes native FSH antigen, Lee Biosolutions 996-11
Buffer	0.9 % NaCl, 0.095 % NaN ₃ as a preservative
IEF Profile	6.6-7.4
Cross Reactivity	FSH α 1 %, FSH β 68 %, LH < 0.2 %, hCG < 0.1 %, TSH < 0.1 %
Specificity	Antibody recognizes human follicle-stimulating hormone (FSH)
Matched Pairs	Anti-h Alpha Subunit 6601 (100066)
Shelf Life	24 months