

Product Datasheet

Human VEGF-A ELISA Kit 650.080

Product Name	Human VEGF-A ELISA Kit
Catalog Number	650.080
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	31.25-1000 pg/ml
Sensitivity	7.9 pg/ml
Tested Applications	ELISA
Brand	Diaclone
Storage	Store at +4°C
Specificity	Specific of human Vascular Endothelial Growth Factor A
Sample Types	Serum, Plasma, Cell culture supernatant