

Product Datasheet

Human CD40L ELISA Kit 650.120

Product Name	Human CD40L ELISA Kit
Catalog Number	650.120
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	0.16-10 ng/ml
Sensitivity	0.062 ng/ml
Tested Applications	ELISA
Alternative Names	TRAP /CD154
Brand	Diaclone
Storage	Store at +4°C
Specificity	Recognizes both natural and recombinant human sCD40L
Sample Types	Serum, Plasma, Cell culture supernatant