

Product Datasheet

Murine IL-1 alpha ELISA Kit 660.010

Product Name	Murine IL-1 alpha ELISA Kit
Catalog Number	660.010
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	7.8-500 pg/ml
Sensitivity	6 pg/ml
Tested Applications	ELISA
Alternative Names	IL-1a
Brand	Diaclone
Storage	Store at +4°C
Specificity	Recognizes both natural and recombinant murine IL-1 α
Sample Types	Serum, Plasma, Cell culture supernatant