

## Product Datasheet

### Murine IL-4 ELISA Kit 660.020

<b>Product Name</b>	Murine IL-4 ELISA Kit
<b>Catalog Number</b>	660.020
<b>Description</b>	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
<b>Range</b>	3.9-250 pg/ml
<b>Sensitivity</b>	3 pg/ml
<b>Tested Applications</b>	ELISA
<b>Brand</b>	Diaclone
<b>Storage</b>	Store at +4°C
<b>Specificity</b>	Recognizes both natural and recombinant murine IL-4
<b>Sample Types</b>	Serum, Plasma, Cell culture supernatant