

Product Datasheet

Murine IL-5 ELISA Kit 660.030

Product Name	Murine IL-5 ELISA Kit
Catalog Number	660.030
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	78-500 pg/ml
Sensitivity	5 pg/ml
Tested Applications	ELISA
Brand	Diaclone
Storage	Store at +4°C
Specificity	Recognizes both natural and recombinant murine IL-5
Sample Types	Serum, Plasma, Cell culture supernatant